

# Ultrasound Inspection Form – Enhanced

## Section 1 – Client Information

Client	
Site	
Date	14/06/2023
WO	JN 1809 ID 12270
Inspector	

## Section 2 – System

System		Inspection Task	Pass/Fail	Notes
Make	GE	Pre-Inspection	Pass	
Model	LOGIQ V2	Physical	Pass	
Serial Number	6032157WX0	Functional	Pass	
Notes	R1.2.1			



## Section 3 – Electrical Safety Test

Safety Tester Reference	A0368	Class & Type	2 BF
Results with Pass/Fail (Pass Limits - See AUE-SOP-018)			
Test	NC	SFC	Pass/Fail
Earth Continuity	N/A	N/A	N/A
Earth Leakage	78 $\mu$ A	N/A	Pass
Enclosure Leakage	13 $\mu$ A	77 $\mu$ A	Pass
Patient Leakage	<4 $\mu$ A	<4 $\mu$ A	Pass
Patient Leakage (F-Type)	N/A	<25 $\mu$ A	Pass
Notes			

## Section 4 – Probe/s

Phantom Reference	A0392
-------------------	-------

Probe ID	11851	Inspection Task	Pass/Fail	Notes
Make	GE	Physical	Pass	
Model	L 6-12 RS	Reverberation & Drop Out	Pass	
Serial Number	734817WX0	Phantom	Pass	

