

Ultrasound Inspection Form – Enhanced

Section 1 – Client Information

Client	
Site	
Date	8/12/22
WO	JN767 ID14696
Inspector	AP

Section 2 – System

System		Inspection Task	Pass/Fail	Notes
Make	Toshiba	Pre-Inspection	Pass	
Model	Aplio 500	Physical	Fail	TGC slider knob missing End Exam button has handwritten label 3 cable roller arm has broken off. DVD drive bezel missing/damaged
Serial Number	WAD1682869	Functional	Pass	
Notes				



Section 3 – Electrical Safety Test

Safety Tester Reference	A0393	Class & Type	I BF
Results with Pass/Fail (Pass Limits - See AUE-SOP-018)			
Test	NC	SFC	Pass/Fail
Earth Continuity	87mΩ	N/A	Pass
Earth Leakage	37μA	N/A	Pass
Enclosure Leakage	<4μA	37μA	Pass
Patient Leakage	<4μA	<4μA	Pass
Patient Leakage (F-Type)	N/A	<25μA	Pass
Notes			

Section 4 – Probe/s

Phantom Reference	A0393
--------------------------	-------

Probe ID	14424	Inspection Task	Pass/Fail	Notes
Make	Toshiba	Physical	Pass	Stock probe used
Model	PLT-1202S	Reverberation & Drop Out	Pass	
Serial Number	99A1683648	Phantom	Pass	

