

Ultrasound Inspection Form – Enhanced

Section 1 – Client Information

Client	
Site	
Date	17/1/2023
WO	JN1031 ID14629
Inspector	

Section 2 – System

System		Inspection Task	Pass/Fail	Notes
Make	GE	Pre-Inspection	Pass	
Model	LOGIQ E9	Physical	Pass	
Serial Number	125761US8	Functional	Fail	Keyboard requires repai.
Notes				

Section 3 – Electrical Safety Test

Safety Tester Reference	A0368	Class & Type	I BF
Results with Pass/Fail (Pass Limits - See AUE-SOP-018)			
Test	NC	SFC	Pass/Fail
Earth Continuity	134mΩ	N/A	Pass
Earth Leakage	206μA	N/A	Pass
Enclosure Leakage	<4μA	204μA	Pass
Patient Leakage	<4μA	<4μA	Pass
Patient Leakage (F-Type)	N/A	<25μA	Pass
Notes			

Section 4 – Probe/s

Phantom Reference	A0393
--------------------------	-------

Probe ID	2567	Inspection Task	Pass/Fail	Notes
Make	GE	Physical	Pass	
Model	IC5-9-D	Reverberation & Drop Out	Pass	
Serial Number	368673WX0	Phantom	Pass	

