

## Probe Report Form – Loan Pool/Stock Probe

<b>Section 1 – Client &amp; Probe Information</b>	
<b>WO &amp; Probe Ref</b>	1-3579 6333
Date	11/12/2018
Probe Make	Philips
Probe Model	X6-1 (explora)
Probe Serial Number	B26889
Probe Part Number	4563561247507
Inspected By	BL

<b>Section 2 – Physical Inspection</b>	
<b>Component</b>	<b>Pass/Fail &amp; Comments</b>
Lens	Pass – Good condition
Probe Casing	Pass – Good condition
Probe Casing Strain Relief	Pass – Good condition
Cable	Pass – Fair condition
Connector Strain Relief	Pass – Good condition
Connector	Pass – Good condition

<b>Section 3 – Electrical Safety</b>		
<b>Test Method</b> (Insert below)	<b>Safety Tester Ref</b>	A0026
Probe Test	<b>Pass/Fail</b>	Pass
Results & Comments (Pass Limit for General & Endocavity Probes <50µA, as per SOP-007)	5.17µA	

<b>Section 4 – Performance</b>		
<b>Test Method</b> (Insert below)	<b>Test Equipment Ref</b>	A0085
System	<b>Pass/Fail</b>	Pass
Comments	Good array performance	

<b>Section 5 – Summary</b>	
Summary	Probe passes inspection & tests
<i>Probe has passed performance and safety assessments and is issued as fit for use, subject to final acceptance assessment upon fit.</i>	







