

# Ultrasound Inspection Form – Enhanced

## Section 1 – Client Information

Client	
Site	
Date	08/03/2024
WO	JN3206 ID18004
Inspector	

## Section 2 – System

System		Inspection Task	Pass/Fail	Notes
Make	Phillips	Pre-Inspection	Pass	
Model	Affiniti CVX	Physical	Pass	
Serial Number	US323R0292	Functional	Pass	
Notes				

## Section 3 – Electrical Safety Test

<b>Safety Tester Reference</b>	A0410	<b>Class &amp; Type</b>	1 BF
<b>Results with Pass/Fail</b> (Pass Limits - See Modality-SOP-018)			
<b>Test</b>	<b>NC</b>	<b>SFC</b>	<b>Pass/Fail</b>
Earth Continuity	268 mΩ	N/A	Pass
Earth Leakage	69 μA	N/A	Pass
Enclosure Leakage	<4μA	55 μA	Pass
Patient Leakage	<4μA	<4μA	Pass
Patient Leakage (F-Type)	N/A	<25μA	Pass
Notes			

## Section 4 – Probe/s

<b>Phantom Reference</b>	A0393
--------------------------	-------

Probe ID	18007	Inspection Task	Pass/Fail	Notes
Make	Phillips	Physical	Pass	
Model	X5-1(CX50)	Reverberation & Drop Out	Pass	
Serial Number	F0085Y	Phantom	Pass	