



Ultrasound Transducer Condition Report

Date Tested 16/10/2020	Carried out By – I. Rutherford
Manufacturer - Philips	Model Number C5-1
Serial Number – B00YXB	
Transducer Type – Convex	
Date of Manufacture - Unknown	
System used to Test – Philips iU22	
Cable – Good	
Lens – Good	
Discolouration - Slight	
Overall Condition – Good	© British Medical Auctions

Airscans –



Ultrasound Transducer Condition Report

Date Tested 16/10/2020	Carried out By – I. Rutherford
Manufacturer - Philips	Model Number – L9-3
Serial Number – 02P455	
Transducer Type – Linjea	
Date of Manufacture - Unknown	
System used to Test – Philips iU22	
Cable – Good	
Lens – Good	
Discolouration –Some	
Overall Condition – Good	

© British Medical Auctions

Airscans –

